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d SDMA IS INCREASED AND CREATININE IS WITHIN THE REFERENCE INTERVAL. IDEXX SDMA Test is a more reliable indicator of kidney function than creatinine because SDMA detects declining kidney function earlier and is not impacted by muscle mass. Creatinine can miss early function loss and be falsely decreased in patients with poor muscle mass. SDMA increases in acute and active injury as well as chronic kidney disease. A complete urinalysis should be performed to evaluate for inappropriate specific gravity, proteinuria and other evidence of kidney disease. If SDMA is 15-19 µg/dL and there is other evidence of kidney disease, action should be taken. If SDMA is 15-19 µg/dL with no other evidence of disease, recheck in 2-4 weeks. If SDMA 120 µg/dL or is persistently 15-19 µg/dL, regardless of whether there is other evidence of kidney disease, action should be taken. For information on recommended actions visit: www.idexx.com/SDMAalgorithm
Note: In puppies, SDMA reference interval is 0-16 µg/dL. SDMA reference interval studies are underway for Kittens and Greyhounds, and results should be interpreted in light of other findings.

Total T4, 11 de julio de 2018

Endocrinology	11.07.18 10:52		
Total T4	24.45	10.3 - 60.49 nmol/L	

b Interpretive ranges:
<0.8 Subnormal
0.8-4.7 Normal
2.3-4.7 Grey zone in old or symptomatic cats
>4.7 Consistent with hyperthyroidism

Cats with subnormal T4 values are almost exclusively euthyroid sick or overtreated for their hyperthyroidism. Older cats with consistent clinical signs and T4 values in the grey zone may have early hyperthyroidism or a concurrent non-thyroidal illness. Hyperthyroidism may be confirmed in these cats by adding on a free T4 or by performing a T3 suppression test. Following treatment with methimazole, T4 values will generally fall within the lower end of the reference range (0.8 - 2.3).

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Urinalysis	14.07.18 15:11	15:06
Collection	Cystocentesis	
Colour	Pale Yellow	
Clarity	Slightly Cloudy	
Specific Gravity	1.022	
pH	5.0	
Urine Protein	30	mg/dL
Glucose	neg	
Ketones	neg	
Blood / Haemoglobin	250	Ery/µL
Bilirubin	neg	
Urobilinogen	norm	
White Blood Cells	1 - 5 /HPF	
Red Blood Cells	6 - 20 /HPF	
Bacteria		
Bacteria, Cocci	Suspect presence*	
Bacteria, Rods	None to rare	
Epithelial Cells		
Squamous Epithelial Cells	None to rare	
Non-Squamous Epithelial Cells	1 - 2 /HPF	
Mucus		
Casts		
Hyaline Casts	None to rare	
Non-Hyaline Casts	None to rare	
Crystals		
Calcium Oxalate Dihydrate Crystals	None to rare	
Struvite Crystals	None to rare	
Unclassified Crystals	None to rare	
Images		

Confirm with one of the following:
Image review
Air-dried, stained cytological preparation ('dry prep')
Urine culture

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
a SDMA IS INCREASED AND CREATININE IS WITHIN THE REFERENCE INTERVAL. IDEXX SDMA Test is a more reliable indicator of kidney function than creatinine because SDMA detects declining kidney function earlier and is not impacted by muscle mass. Creatinine can miss early function loss and be falsely decreased in patients with poor muscle mass. SDMA increases in acute and active injury as well as chronic kidney disease. A complete urinalysis should be performed to evaluate for inappropriate specific gravity, proteinuria and other evidence of kidney disease. If SDMA is 15-19 µg/dL and there is other evidence of kidney disease, action should be taken. If SDMA is 15-19 µg/dL with no other evidence of disease, recheck in 2-4 weeks. If SDMA ≥20 µg/dL or is persistently 15-19 µg/dL, regardless of whether there is other evidence of kidney disease, action should be taken. For information on recommended actions visit: www.idexx.com/SDMAAlgorithm
Note: In puppies, SDMA reference interval is 0-16 µg/dL. SDMA reference interval studies are underway for Kittens and Greyhounds, and results should be interpreted in light of other findings.

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